



The threads are classified in ABSORBABLE and NON ABSORBABLE. They can be from natural and synthetic origin.



> **NEEDLE WITH JAPANESE STEEL AISI 300 SERIES**

- ☑ Silicone, which allows a smooth passage through tissues.
- ☑ Designed to keep the edge after several passes.
- ☑ Needle of suitable size for each caliber.


> **HIGH QUALITY THREAD**

- ☑ Uniform and homogeneous strand, which prevents tissue damage.
- ☑ High tensile strength.
- ☑ Flexible, facilitates handling and safety in performing the knots.
- ☑ Thread with proper length for each type of intervention.

> **BENEFITS FOR THE PATIENT**

- ☑ Less tissue reaction and smooth passage through tissue due to the excellent overall silicone edge and needles.
- ☑ Perfect attachment between the needle and the thread generates less trauma in the patient.

> **GRAPHIC OF OUR NEEDLES**



> **NEEDLE POINTS**

1/2 CR 1/2 Circle tip 3/8 CR 3/8 Circle tip
 1/2 CC 1/2 Cutting circle tip 3/8 CC 3/8 Cutting circle tip

Type	Point	Shape	Pattern
Cutting Point			
Taper Point			
Spatula Point			

ABSORBABLE SUTURES

Natural

Plain Gut

2-0
USP
(3.5 Metric)
28" (70 cm)
MR 15

General Surgery / Cirugía General
PLAIN GUT
CATGUT SIMPLE
Absorbable / Surgical Suture
Absorbible / Sutura Quirúrgica
(STERILE) EO

1/2 ∅ 15 mm
Taper Point / Punta Redonda
R. San.: DM026N

TYPE OF ABSORPTION:
By phagocytosis

STRUCTURE:
Monofilament

COLOR:
Cream, bright yellow

USP SIZE:
6-0 to 3.

Chromic Gut

2-0
USP
(3.5 Metric)
28" (70 cm)
MR 15

General Surgery / Cirugía General
CHROMIC GUT
CATGUT CRÓMICO
Absorbable / Surgical Suture
Absorbible / Sutura Quirúrgica
(STERILE) EO

1/2 ∅ 15 mm
Taper Point / Punta Redonda
R. San.: DM026N

TYPE OF ABSORPTION:
By phagocytosis

STRUCTURE:
Monofilament

COLOR:
Bright brown, dark brown

USP SIZE:
6-0 to 3.

Synthetic

Polyglycolic Acid

2-0
USP
(3 Metric)
28" (70 cm)
MR 15

General Surgery / Cirugía General
POLYGLYCOLIC ACID
ÁCIDO POLIGLICÓLICO
Absorbable / Surgical Suture
Absorbible / Sutura Quirúrgica
(STERILE) EO

1/2 ∅ 15 mm
Taper Point / Punta Redonda
R. San.: DM026N

TYPE OF ABSORPTION:
By hydrolysis

STRUCTURE:
Braided
Multifilament

COLOR:
Violet

USP SIZE:
8-0 to 2.

Polyglactin

2-0
USP
(3 Metric)
28" (70 cm)
MR 15

General Surgery / Cirugía General
POLYGLACTIN
POLYGLACTIN
Absorbable / Surgical Suture
Absorbible / Sutura Quirúrgica
(STERILE) EO

1/2 ∅ 15 mm
Taper Point / Punta Redonda
R. San.: DM026N

TYPE OF ABSORPTION:
By hydrolysis

STRUCTURE:
Braided
Multifilament

COLOR:
Violet

USP SIZE:
8-0 to 2.

Polydioxanone

2-0
USP
(3 Metric)
28" (70 cm)
MR 15

General Surgery / Cirugía General
POLYDIOXANONE
POLIDIOXANONA
Absorbable / Surgical Suture
Absorbible / Sutura Quirúrgica
(STERILE) EO

1/2 ∅ 15 mm
Taper Point / Punta Redonda
R. San.: DM026N

TYPE OF ABSORPTION:
By hydrolysis

STRUCTURE:
Monofilament

COLOR:
Purple

USP SIZE:
7-0 to 2

Poliglecaprone (Unicryl)

1
USP
(4 Metric)
28" (70 cm)
MR 40

General Surgery / Cirugía General
POLYGLECAPRONE
POLIGLECAPRONE
Absorbable / Surgical Suture
Absorbible / Sutura Quirúrgica
(STERILE) EO

1/2 ∅ 40 mm
Taper Point / Punta redonda
R. San.:

TYPE OF ABSORPTION:
By hydrolysis

STRUCTURE:
Monofilament

COLOR:
Violet, undyed

USP SIZE:
5-0 to 1

NON ABSORBABLES SUTURES

Synthetic

Polypropylene

3-0
USP
(2 Metric)
28" (70 cm)
MR 25

General Surgery / Cirugía General
POLYPROPYLENE
POLIPROPILENO
Non Absorbable / Surgical Suture
No Absorbible / Sutura Quirúrgica
(STERILE) EO

1/2 ∅ 25 mm
Taper Point / Punta redonda
R. San.:

TYPE OF ABSORPTION:
Non-absorbable

STRUCTURE:
Monofilament

COLOR:
Blue

USP SIZE:
10-0 to 1.

Polyester

2-0
USP
(3 Metric)
28" (70 cm)
MR 25

General Surgery / Cirugía General
POLYESTER
POLIESTER
Non Absorbable / Surgical Suture
No Absorbible / Sutura Quirúrgica
(STERILE) EO

1/2 ∅ 25 mm
Taper Point / Punta redonda
R. San.:

TYPE OF ABSORPTION:
Non-absorbable

STRUCTURE:
Braided
Multifilament

COLOR:
Green

USP SIZE:
6-0 a 5.

Virgin Silk

8-0
USP
(0.4 METRIC)
20PSS3
6.4 mm

General Surgery / Cirugía General
VIRGIN SILK
SEDA VIRGEN
Non Absorbable / Surgical Suture
No Absorbible / Sutura Quirúrgica
(STERILE) EO

3/8 Circle Spatulated / Círculo Espalada
R. San.:

TYPE OF ABSORPTION:
Non-absorbable

STRUCTURE:
Twisted
Multifilament

COLOR:
Blue, black

USP SIZE:
8-0 to 7-0.

Black Braided Silk

3-0
USP
(2 Metric)
28" (70 cm)
TC 20

General Surgery / Cirugía General
BLACK BRAIDED SILK
SEDA NEGRA TRENZADA
Non Absorbable / Surgical Suture
No Absorbible / Sutura Quirúrgica
(STERILE) EO

3/8 Reverse Cutting / Reverso Cortante
R. San.:

TYPE OF ABSORPTION:
Non-absorbable

STRUCTURE:
Braided
Multifilament

COLOR:
Black, White

USP SIZE:
7-0 to 2.

Black Nylon

2-0
USP
(3 Metric)
28" (70 cm)
TC 25

General Surgery / Cirugía General
MONOFILAMENT NYLON
NYLON MONOFILAMENTO
Non Absorbable / Surgical Suture
No Absorbible / Sutura Quirúrgica
(STERILE) EO

3/8 Reverse Cutting / Reverso Cortante
R. San.:

TYPE OF ABSORPTION:
Non-absorbable

STRUCTURE:
Monofilament

COLOR:
Black

USP SIZE:
10-0 to 2.

Recommendations:

Open only sutures used for the procedure. With gentle traction to straighten the sutures. Do not tighten, not rub, not abruptly pull the needles. Do not press or squeeze the threads of sutures with surgical instruments.

Warnings:

Do not resterilize. Discard open unused sutures. Store at a temperature not exceeding of 27°C.